

Checking Your Blood Pressure

Before You Begin

- ☐ **No coffee**
in the last 30 minutes.



- ☐ **No tobacco**
in the last 30 minutes.



- ☐ **No meals or exercise**
within 2 hours.




- ☐ **Visit the rest room**
A full bladder can affect the reading.






Checking Your Blood Pressure

Getting Ready




- ◆ Insert the AC adapter plug into the AC adapter jack (1) and plug into an outlet (2).
- ◆ Select the correct blood pressure cuff size. Wipe the inside of the cuff with a disinfectant wipe.
- ◆ Insert the air plug securely into the main unit (3).

- ◆ Sit quietly in a chair for 5 minutes before measuring.
- ◆ Keep your back supported and feet flat on the ground.
- ◆ Use the same arm each time you take your blood pressure.
- ◆ Remove tight-fitting clothing from your upper arm.
- ◆ The bottom of the cuff should be about a half inch above your elbow.




Checking Your Blood Pressure

Taking Your Blood Pressure



- ◆ Support your arm at the level of your heart on a table or high armrest.
- ◆ Stay in a seated position.
- ◆ Avoid talking when taking a measurement.
- ◆ **Take 3 readings, each 2 minutes apart.**
- ◆ Press the START/STOP button on the monitor. Display symbols will appear on the screen.



- ◆ The cuff will inflate on its own and feel tight for a few seconds.
- ◆ The cuff will relax and display the reading. Record the reading.
- ◆ Press STOP.
- ◆ Wait two minutes. Repeat these steps for a total of three readings.

Checking Your Blood Pressure

What The Numbers Mean

Know your numbers – write them down!

- ◆ Write down blood pressure and pulse rate numbers displayed on the screen.
- ◆ **If the reading is very different from what you expected, please see your doctor to confirm the reading.**
- ◆ Also, make a note on the comment sheet so that we can check out any problems with the equipment.

Within Healthy Range: Less than 120 <i>systolic</i> Less than 80 <i>diastolic</i>	Excellent! This is right where it should be and is a normal blood pressure.
Approaching High Risk: 120–139 <i>systolic</i> 80–89 <i>diastolic</i>	This is getting up there. You should recheck it within 2 weeks and begin lifestyle changes that can help keep it from becoming full blown hypertension.
In High Range: See a doctor soon 140–170 <i>systolic</i> 90–100 <i>diastolic</i>	Your blood pressure is high. You should see a doctor.
High Reading: Seek medical care now 170 or higher <i>systolic</i> 100 or higher <i>diastolic</i>	Your blood pressure is very high. Seek medical care now!

Created by the Heart Disease & Stroke Prevention Program, Washington State Department of Health



Dear Colleague,

The Washington State Department of Health (DOH) provides print-ready files (PDFs) of health education materials. To ensure that the original quality of the piece is maintained, please read and follow the instructions below and the specifications included for professional printing.

- **Use the latest version.** DOH materials are developed using the most current information available, are checked for clinical accuracy, and are field tested with the intended audience to ensure they are clear and readable. DOH programs make periodic revisions to educational materials, so please check this web site to be sure you have the latest version. DOH assumes no responsibility for the use of this material or for any errors or omissions.
- **Do not alter.** We are providing this artwork with the understanding that it will be printed without alterations and copies will be free to the public. Do not edit the text or use illustrations or photographs for other purposes without first contacting us. Please do not alter or remove the DOH logo, publication number or revision date. If you want to use a part of this publication for other purposes, contact the Office of Health Promotion first.
- **For quality reproduction:** Low resolution PDF files are intended for black and white or color desktop printers. They work best if you are making only one or two copies. High resolution PDF files are intended for reproducing large quantities and are set up for use by professional offset print shops. The high resolution files also include detailed printing specifications. Please match them as closely as possible and insist on the best possible quality for all reproductions.

If you have questions, contact:
Office of Health Promotion
P.O. Box 47833 Olympia, WA 98504-7833
(360) 236-3736

Sincerely,
Health Education Resource Exchange Web Team

P R I N T I N G S P E C I F I C A T I O N S

Title: **Checking Your Blood Pressure** (4 posters)

Size: 11 x 17

Paper stock: 80# Mohawk 50/10 Matte Cover, white. Matte protective coating 1 side.

Ink color: 4-color process

DOH Pub #: 345-267